

AutoScratch™

Wound Making Tool

GENERAL	
Microplate types	96- and 24-well microplates (BioTek PNs 02638 and 02635)
Automated functions	Wound making Cleaning process
Scratch pins and manifolds	User removable, replaceable pins User-exchangeable manifolds
Wound Parameters	Entire wound visible within single 4x image taken on Cytation or Lionheart
Cleaning / sterilization	Built-in automated scratch pin cleaning procedure Built-in cleaning reagent troughs removable for chemical sterilization
PHYSICAL CHARACTERISTICS	
Connectivity	1 USB 2.0 port for software download
Dimensions	12" L x 7.5" W x 6.75" H (30.4 cm x 19.1 cm x 17.2 cm)
Weight	8 lbs (3.6 kg)
Power consumption	40 W maximum
REGULATORY	
Regulatory	CE and TUV marked. RoHS compliant.

Specifications are subject to change. Performance values represent the average observed factory test values.