

Imaging System Controller (PN 1760012)

Computer for Imaging Applications

COMPUTER CONFIGURATION	
HP Part Description	HP ProDesk 400 G3 Desktop Mini Business PC
Operating System	MS Windows 10 Pro, 64-Bit
Processor	Intel 7500T 7GEN i5 – 4 Core (6MB cache, 4 threads)
Memory	8 GB DDR42400 (1x8 GB)
Storage	512GB PCIe NVM3 TLC SSD
Input Device	USB Business slim keyboard and mouse
Power	65W DM Ext power adapter (US adapter)
Warranty	3/3/3 DM Warranty
Ports & Connectors	(4) USB 3.1, (2) USB 2.0; - (Keyboard and Mouse), Serial Port, HMDI, RJ-45 Network, Dual-Mode DisplayPort
Restore Media	One bootable USB Thumb drive containing HP recovery media.
MONITOR (BTI PN 02346)	
Display size (diagonal)	23.8" (60.45 cm)
Native resolution	1920 x 1080 @60 Hz (FHD)
Input connectors	1 VGA, 1 DisplayPort 1.2 (with HDCP support), 1 HDMI (with HDCP support)
Weight	9.47 lb (4.3 kg) <i>(With stand)</i>
Dimensions (W x D x H)	With stand: 22.49 x 7.53 x 16.39 in (57.12 x 19.12 x 41.64 cm) Without stand: 22.49 x 2.17 x 14 in (57.12 x 5.5 x 35.55 cm)

Rev. IMSPC040418